


## Bio-Decontamination with H<sub>2</sub>O<sub>2</sub> – Vapor in Pharmaceutical Industry

Dr. Johannes Rauschnabel,  
Robert Bosch GmbH



SBM Expert's Congress 13<sup>th</sup>-14<sup>th</sup> of September, 2007

Packaging Technology

1


PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

## Structure

- Introduction
- “VHP” for Isolators
- Results
- More Applications
- Summary



Packaging Technology

2


PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

### Structure

- Introduction
  - Chemistry / Biocidal Effect
  - Challenges
  - Regulatory
- “VHP” for Isolators
- Results
- More Applications
- Summary



### Introduction



- Hydrogen Peroxide is an old sterilant – well established in food industry, where liquid H<sub>2</sub>O<sub>2</sub> is applied for disinfection of packaging material and plastic bottles.
- H<sub>2</sub>O<sub>2</sub> has a broad antimicrobial activity spectrum and shows low toxicity.
- Residues of hydrogen peroxide decompose into neutral species, which underlines the environmental friendliness of this agent:
$$2 \text{H}_2\text{O}_2 \rightleftharpoons 2 \text{H}_2\text{O} + \text{O}_2$$
- But: hydrogen peroxide can decompose with fast kinetics – under some circumstances it can explode at high concentrations.
- For complex geometries disinfection with liquid peroxide has limits. For this case H<sub>2</sub>O<sub>2</sub> in the vapor state is widely used.




## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Chemistry / Biocidal Effect

- The hydrogen peroxide molecule has a unique geometry with a high energy angle, which makes splitting of O-O-bond easy:



- The decomposition can be catalyzed by metallic surfaces (Cu, Pb, Mg...) – which can be used to decompose after use. On the other side this makes the need for noble surfaces in the equipment to be disinfected.



Packaging Technology

5

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Chemistry / Biocidal Effect

- The biocidal effect of hydrogen peroxide comes from the oxygen atom/radical, which is released under influence of the enzyme catalase at the microbotic surface. The oxidative power of oxygen radicals cause irreversible damage to enzymic systems and DNA.
- A minimum of humidity at the surface to be treated helps to promote this metabolismism – especially bacterial spores need a basic humidity to show any biochemical activity.



Packaging Technology

6

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Structure

- Introduction
  - Chemistry / Biocidal Effect
  - Challenges
  - Regulatory
- “VHP” for Isolators
- Results
- More Applications
- Summary



Packaging Technology

7

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Challenges

- Vaporized Hydrogen Peroxide is used to disinfect bottles and pharmaceutical equipment interior, such as isolators (mini-cleanrooms on top of machinery).
- Pioneer in that technique was AMSCO, later bought by Steris Co., with their VHP® technology back in the 80ies.
- Disadvantage of vaporizing H<sub>2</sub>O<sub>2</sub> are the two physical states present in the system, when using diluted H<sub>2</sub>O<sub>2</sub> for evaporation. This results in different concentrations and thereby in different inactivation efficacy.



Packaging Technology

8


PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**


### Challenges

- Slow evaporation of e.g. 35% H<sub>2</sub>O<sub>2</sub> separates the both components (“distillation”) during transition into vapor. Condensation on the surfaces makes the next separation (“concentration change”).
- Sterilisation means absence of viable microbotics. In complex geometries this is a challenge, because homogeneous inactivation is required. This can be questioned with vapor/liquid mixtures.
- H<sub>2</sub>O<sub>2</sub> vapor is known as *bio-decontaminating, not sterilizing* surfaces – though 6 log inactivation can be achieved/proven.



### Structure

- **Introduction**
  - Chemistry / Biocidal Effect
  - Challenges
  - **Regulatory**
- “VHP” for Isolators
- Results
- More Applications
- Summary





### Regulatory

- In “Guidance for Industry: Sterile Drug Products Produced by Aseptic Processing — Current Good Manufacturing Practice” published by FDA in September 2004, within Appendix 1: Aseptic Processing Isolators the need for decontaminated surfaces is explained:



#### *“1. Surface Exposure*

*Decontamination procedures should ensure full exposure of all isolator surfaces to the chemical agent. The capability of a decontaminant to penetrate obstructed or covered surfaces is limited....”*



### Regulatory

- Continuation “Aseptic Guidance”:

#### *“2. Efficacy*

*The decontamination method should render the inner surfaces of the isolator free of viable microorganisms. Multiple available vaporized agents are suitable for achieving decontamination.... Cycles should be developed with an appropriate margin of extra kill to provide confidence in robustness of the decontamination processes. Normally, a four- to six-log reduction can be justified depending on the application.”*

#### *“E. Filling Line Sterilization*

*... Where decontamination methods are used to render certain product contact surfaces free of viable organisms, a minimum of a six-log reduction should be demonstrated using a suitable biological indicator.”*



### Regulatory

- Continuation “Aseptic Guidance”:  
“E. Filling Line Sterilization  
... Where decontamination methods are used to render certain product contact surfaces free of viable organisms, a minimum of a six-log reduction should be demonstrated using a suitable biological indicator.”
- In EC GUIDE to Good Manufacturing Practice, Annex 1 a list of sterilization techniques (heat, moist heat, radiation, ethylene oxide...) does not include hydrogen peroxide. The description of sterilization/decontamination of isolators does not specifically discuss issues coming from bio-decontamination instead sterilization.



### Regulatory

- In EC GUIDE to Good Manufacturing Practice, Annex 1 a list of sterilization techniques (heat, moist heat, radiation, ethylene oxide...) does not include hydrogen peroxide. The description of sterilization/decontamination of isolators does not specifically discuss issues coming from bio-decontamination instead “sterilization”.



### Regulatory

- PIC/S Recommendation “Isolators used for Aseptic Processing and Sterility Testing”; Chapter “An expansion of the detailed points to be considered for the implementation of the principles to isolators subjected to a sporicidal process” collects 14 points for the sporicidal process. Point 9.4.1 mentions H<sub>2</sub>O<sub>2</sub> bio-decontamination:


*“... The agent selected for gas generation should be sporicidal. The agent used for gas generation or other means of application should be capable of rapidly killing bacterial endospores, fungal spores and vegetative microorganisms.*

*... Peracetic acid, hydrogen peroxide and formaldehyde are used. ...”*



### Structure

- Introduction
- **“VHP” for Isolators**
- Results
- More Applications
- Summary





## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

- Vaporized hydrogen peroxide is widely used to decontaminate the inside of pharmaceutical isolators. The good old pioneering days applied a mixture of steam and H<sub>2</sub>O<sub>2</sub> to be on the safe side, but validation was a nightmare.
- Nowadays several H<sub>2</sub>O<sub>2</sub> gas generators are on the market – different philosophy, but more or less reliable technology.
- Cycle development and validation has become a routine job.



Packaging Technology

17

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

- Isolator design principle



Packaging Technology

18

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

Example (isolated vial/bottle fill line in Austria)



Packaging Technology

19

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators


#### → Closed loop system:

Liquid H<sub>2</sub>O<sub>2</sub> is dropped on a heated plate inside the filter level of an isolator. Recirculation and continuous injection of peroxide goes to saturation of atmosphere including strong visible condensation



#### → Open loop system:

Liquid hydrogen peroxide is evaporated outside of the isolator in a gas generator; which carries the vapor with air into the isolator, where it is recirculated. Some o-l-systems perform “flash evaporation”.



Packaging Technology

20

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators


#### → Closed loop system:

- Simple air management
- Vapor passes filters first
- Incomplete evaporation
- Limited injection rate
- Saturation/condensation soon



#### → Open loop system:

- Complete (“flash”) evaporation
- Controlled atmosphere
- Capability for high gas conc.
- Air management more complex
- More sophisticated generator



Packaging Technology

21

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

#### Preparation (Phase I)

- Expanding gloves, load monitoring goods
- Preparation of isolator atmosphere
- Preheating of peroxide duct work

#### Conditioning (Phase II)


- Injection of H<sub>2</sub>O<sub>2</sub> with a high rate to achieve a high gas concentration

#### Biodecontamination (Phase III)

- Injection with a lower rate to replace inactivated peroxide and keep level

#### Aeration (Phase IV)

- Removal of H<sub>2</sub>O<sub>2</sub> from injection system
- Removal of peroxide from all isolator surfaces including filters.



Packaging Technology

22

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**


## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

→ Bosch SafeVAP (open loop system)

Features:

- Venturi transport
- 3-step evaporation
- Inj. rates up to 80g / min
- Controlled atmosphere



**SafeVAP = Safe Vacuum Assisted Peroxide Evaporation**

Packaging Technology

23

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

→ Description



Packaging Technology

24

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

#### → Description



Packaging Technology

25

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### “VHP” for Isolators

- Cycle Development for H<sub>2</sub>O<sub>2</sub> bio-decontamination include
  - OQ-studies
  - tests with chemical indicators (CI's)
  - Tests with bio indicators (BI's), which are stainless steel coupons inoculated with 6 log of the most resistant viable (spores of *Geobacillus stearothermophilus*).
- To ease handling (and to put an extra challenge on the biological test system), BI's wrapped in Tyvek became industrial standard.



Packaging Technology


26

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

### “VHP” for Isolators



### Structure

- Introduction
- “VHP” for Isolators
- **Results**
  - **Recent Studies**
  - Results
- More Applications
- Summary





## Recent Studies

- Many studies were performed at Bosch during the last 3 years (PhD thesis). Basis for these studies was a reference isolator equipped with the SafeVAP system and placed in the Bosch PharmaLab.



Packaging Technology

29



PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

## Recent Studies

- Material of Construction Study:  
Goal of this study was to learn the influence of different surfaces (material, finish) on inactivation efficacy of bioburden.
- Humidity Study  
Influence of relative humidity at different H<sub>2</sub>O<sub>2</sub> gas concentration levels was studied and compared with dew point analysis.



Packaging Technology

30

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Structure

- Introduction
- “VHP” for Isolators
- **Results**
  - Recent Studies
  - **Results**
- More Applications
- Summary



Packaging Technology

31

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

Objective of material study:

- Examine influence of different construction material surfaces on inactivation of spores
- Determine variable kill kinetics (D-Value)
- Test organism: spores of *Geobacillus stearothermophilus* ATCC 12980 with population of > 1\*10<sup>6</sup> spores per BI.



Exposure time [min]	0	5	10	15	20	25
Results	+	+	+	-	-	-
	+	+	-	-	-	-
	+	+	-	-	-	-
	Survival		Fractional	Killing		

Packaging Technology

32

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

Results on some materials (50 different surfaces tested)



Packaging Technology

33

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

Results on standard control (commercial Bio indicator made from SSL)



Packaging Technology

34

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

### Results

Summary of the study:

- Most materials tested show good inactivation properties.
- Plastic/elastomeric materials have generally higher D-values compared with stainless steel and glass.
- Standard commercial available BI's show a certain spread of D-values (0,2 min to 0,9 min – late positives included in the calculation), which represents microbiological accuracy.
- A correlation between wetting properties, surface roughness and the entire D-values can be recognized.
- Critical materials of construction should either be replaced or need prolonging of the cycle time.



### Results

The results were recently published:

#### **Suitability of Different Construction Materials for Use in Aseptic Processing Environments Decontaminated with Gaseous Hydrogen Peroxide**

BEATRIZ UNGER,<sup>1</sup> UTA RAUSCHNABEL,<sup>1</sup> BERTHOLD DÜTHORN,<sup>1</sup> VOLKER KOTTKE,<sup>2</sup> CHRISTIAN HERTEL,<sup>2</sup> AND JOHANNES RAUSCHNABEL<sup>1</sup>

<sup>1</sup>Robert Bosch GmbH, Blaufelder Str. 45, 74564 Crailsheim, Germany, Phone: +49 7951 402-471, Fax: +49 711 811-5180358, E-mail: [beatriz.unger@de.bosch.com](mailto:beatriz.unger@de.bosch.com); <sup>2</sup>University of Hohenheim, 70599 Stuttgart, Germany

*PDA journal, Vol. 61, No. 4, July-August 2007*




## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

Procedure for Humidity study:

- Separate evaporation of water and hydrogen peroxide (35%)
- Control of humidity level with NIR probe, dew point sensor and standard humidity sensor
- D - value determination of commercial BI's (SSL)
- Investigation on inactivation in hindered surfaces with high aspect ratio



Packaging Technology

37

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

Details:



Packaging Technology

38

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

400 ppm:



Packaging Technology

39

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

600 ppm



Packaging Technology

40

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**



## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

800 ppm:



Packaging Technology

41

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.




**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Results

Comparison  
 Wrapped vs.  
 Nude BI's :



Packaging Technology


42

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

### Results



### Results

#### Summary:


- The higher the cgas phase concentration, the more independent is the inactivation effect from the humidity level.
- At lower concentrations the same kill as at higher concentration can be achieved with higher humidity (condensation)
- The faster the concentration rises during conditioning phase (flash evaporation), the less peroxide is needed to achieve a good kill.
- (Nano/micro) condensation occurs at all humidity levels. Visible concentration is not needed to get short inactivation times.



## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### Structure

- Introduction
- “VHP” for Isolators
- Results
- **More Applications**
- Summary



Packaging Technology

45

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor

### More Applications

- Vaporized H<sub>2</sub>O<sub>2</sub> is applied with all kinds of isolators: fill/finish isolators, sterility testing isolators, hospital isolators...
- A new trend is fumigation of cleanrooms and animal testing facilities with hydrogen peroxide vapor.
- H<sub>2</sub>O<sub>2</sub> in vapor state could potentially replace alcoholic spray&pray procedures in material transfer chambers and material locks for thermolabile products.
- Bio-decontamination of tubs with presterilized, nested syringes can also be performed with H<sub>2</sub>O<sub>2</sub>.



Packaging Technology

46


PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**

### Structure

- Introduction
- “VHP” for Isolators
- Results
- More Applications
- **Summary**



### Summary

- Hydrogen peroxide vapor is widely used in pharmaceutical industry. One major application is disinfection of isolator interior.
- Closed and open loop systems are established in the market, which are technically robust and validatable.
- Inactivation effect on bioburden differs depending on the material of construction (material and finish).
- Gas concentration and humidity effect dominate the kill efficacy: high gas concentrations through flash evaporation help to achieve efficient processes.
- New applications of H<sub>2</sub>O<sub>2</sub> vapor include cleanroom fumigation and decontamination of thermolabile materials in transfer chambers,



## Bio-decontamination with H<sub>2</sub>O<sub>2</sub> - Vapor



*johannes.rauschnabel@de.bosch.com*

Packaging Technology

49

PA-PH/EPP | 10/09/2007 | © Robert Bosch GmbH 2007. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



**BOSCH**